

MED4ELITE



Description:

Game Ready introduces a breakthrough, multi-modality, dual-user recovery system. The new touchscreen-controlled Med4 Elite™ integrates iceless cold, controllable heat, rapid contrast, and active compression therapies like never before. Now, physical therapists, athletic trainers, physicians, and diverse users have comprehensive, flexible, and clinically proven treatment options for diverse injuries, post-surgery recovery, and different stages of rehabilitation. All to help get people back to being their best. This is recovery re-engineered.

TECHNICAL SPECS

- **Dimensions:** 32.5" L x 24.75" W x 43" H (83cm L x 63cm W x 109cm H)
- **Weight:** [empty] 172 lbs (78kg)
- **Hose Connections:** 2 ports, 2 holsters
- **Wheels:** 4 (lockable)

Four Therapeutic Modalities, One Revolutionary System

The Med4 Elite multi-modality therapy unit can improve outcomes effectively and safely with the most innovative and advanced recovery technology available, while streamlining professional workflow by enabling the simultaneous treatment of one or two patients or multiple injuries - all with targeted efficiency.

The system integrates cryotherapy, thermotherapy, rapid contrast therapy (RCT), and intermittent pneumatic compression therapy (IPC) in one device. A touchscreen makes control easy. And four lockable wheels facilitate optimal positioning.

HEAT

- Set temperature UP TO 113°F/45°C.
- Creates its own Heat Therapy
- Set total Therapy time.
- Set Compression level to Low or None.

COLD

- Set temperature, Down To 38°F/3°C.
- Set total Therapy time.
- Set Compression level to Low, Medium-Low, Medium, High, or None.
- Snooze feature, up to 6 cycles.

RAPID CONTRAST

- capable of a hot-cold changeover rate of less than one minute for clean, rapid contrast delivery
- Set Heat and Cold Cycle time and total Therapy time.
- Set Compression level for Cold to Low, Medium-Low, Medium, High, or None.
- Start Rapid Contrast Therapy with either Heat or Cold Therapy.

COMPRESSION-ONLY

- Set total Therapy time.
- Set Compression level to Low, Medium-Low, Medium, or High.

Supported by:

Elite Equipments

✉ support@eliteequipments.com

🌐 <https://eliteequipments.com>

Product Support:

For installation, servicing, or technical assistance, please contact our support team or visit our website.